



HAAS SOFTWARE

Mac OS Software

WWW.HAAS-SOFTWARE.COM

T (888)245-5214


RadioVision User Manual

Haas Software

GETTING STARTED

• Installation

1. Open the RadioVision installer, RadioVision.pkg
2. If you do not have PostgreSQL installed in the /Library/PostgreSQL folder, choose to install PostgreSQL. If you do have PostgreSQL installed already, this option will be disabled.
3. Choosing to install the demo database will create a file called RadioVision.back in your applications folder, which can be imported using the “Restore From Backup...” command in the file menu (see the section on the File Menu at the end of the manual).
4. RadioVision will be installed in your Applications folder.



For new installations, all four packages should be installed. RadioVision will not run without PostgreSQL and the driver support files installed on your computer.


GETTING STARTED

• Opening the RadioVision application.

1. SINGLE CLICK on the RadioVision icon on the dock or DOUBLE CLICK on the RadioVision icon in the applications folder.




2. A registration message will appear when you open RadioVision if you have not licensed the program. You will have a 15 day grace period to complete the licensing process. If the grace is exhausted, RadioVision will default to the demo mode in which image saving is disallowed. If you want to bypass the licensing process click the continue button and proceed to #5.




** If you do not have an email account configured on your server, choose "Manual Request" and call your technical support with the displayed hardware ID; your tech support will provide a serial number that you will enter manually.

3. Click on the Request License button and a new mail window will open that is addressed to Haas Software. The email will contain a file. Be sure that your email address is in the “from” field.
4. You will receive an email with the license file.
 - 4.1. On the server, open the email and double click on the license file. RadioVision will open and the licensing process will be completed. No further action will be required for the client computers.



5. Radiovision is the default username an administrative account, you should change the password for this account and you can add other accounts at this time.
LOG IN NAMES AND PASSWORDS MUST BE LOWER CASE LETTERS OR NUMBERS ONLY !
6. If this is the first time you opened RadioVision on your network and no server has been established, the program will begin searching for the server and you will need to click on the **Use Local** button. **No Active Server** will be indicated in the lower left hand corner of the application window. Selecting the **Use Local** button establishes this computer as the RadioVision server for your network. When installing RadioVision on subsequent computers, the server will automatically be located.

NOTE: THE COMPUTER THAT IS ACTING AS THE SERVER MUST BE ON AND RUNNING RADIOVISION FOR CLIENT COMPUTERS TO LOCATE THE SERVER.




• Setting up user accounts (OPTIONAL)

Only an administrator can add a user account.

An administrator has delete privileges.



1. Go to the Menu Bar and choose [Users](#).
2. A drop down list will appear and you will click on [Manage](#).



3. [radiovision](#) will be one of the usernames and can not be deleted. The password for the radiovision username is radiovision (we recommend you change this password).
4. Click on the [Update](#) button.
5. To add a user, simply highlight the new user and type in the new username and password. The new user is added to the database when you click on the update button.

NOTE: AT LEAST ONE USER NEEDS TO HAVE ADMIN PRIVILEGES


NOTE: RADIOVISION AND THE NAME OF THE USER ACCOUNT USED TO INSTALL RADIOVISION ARE TWO USERS THAT ARE AUTOMATICALLY CREATED

A screenshot of the 'Manage Users' dialog box. The 'User Name' list on the left has 'radiovision' selected. The form fields on the right are: Username: 'radiovision', Password: '.....', New Password: '....', Re-type password: '....', and 'Has Admin Privilege' is checked. The 'Update' button is highlighted with a mouse cursor.A screenshot of the 'Manage Users' dialog box. The 'User Name' list on the left has 'New User' selected. The form fields on the right are: Username: 'John', Password: '.....', New Password: '....', Re-type password: '....', and 'Has Admin Privilege' is checked. The 'Update' button is highlighted with a mouse cursor.


• Adding patients to the database.

After you have logged into RadioVision you will be able to add patients to your database.

1. Single click on the **Select** button.
2. A sheet will drop down from the tool bar and you will need to single click on the **Add New** button.



3. Enter the name of the patient
4. Filling in the Birthdate and Social Security fields is optional for more accurate data collection.
5. Single click on the **Save** button.



6. The name of the new patient will be highlighted in the list.
7. To open the new patient's chart, **DOUBLE CLICK** on their name in the list.

•RadioVision Diagram

Menu Bar

Tool Bars


Detail View

Expanded View

History View

Dock

RadioVision Icon




• Taking X-Rays


Once you are in a patient's chart, you are ready to begin taking x-rays

1. To take an x-ray, single click on the **Take Xray** icon on the tool bar.


NOTE: THE PATIENT'S NAME YOU HAVE SELECTED SHOULD APPEAR ON THE TOOLBAR



- 1.1. The sensor needs to be plugged into the computer at this time. If no sensor is plugged in, an error will appear. Plug in the sensor and click on the **OK** button.




2. A sheet will drop down from the tool bar and you will need to choose the type of x-ray you would like to take. Single click on the **layout type**
3. A sheet will drop down you can select the layout you wish to take. The date is set to the current date and is not editable.
 - 3.1. If you are taking a single film, choose the correct orientation (vertical or horizontal). Then, select the teeth from the diagram that will be included in the film. You may select more than one tooth.
 - 3.2. If you are taking a layout series there is no need to select the teeth that will be included.
4. Ensure that your sensor hardware is selected in the "sensor" list on the upper right.
5. Once you have selected the layout, single click on the **OK** button.




• Taking a single X-Ray

Once you have selected the orientation, teeth and **OK** button, the exposure screen will appear.


1. There will be a **red** outline around the x-ray placeholder.
2. The outline will change to **green** and you will hear a “ding” when the sensor is active and ready to be exposed.



3. Position your sensor.
4. Align X-Ray tube.
5. Push X-Ray button. You should hear another “ding” indicating that the exposure was detected.
6. Your image should appear on your screen within 2 to 6 seconds, depending on the particular sensor you’re using



7. Save the image by single clicking on the **Save** icon on the tool bar.
8. After the images have been saved, they will appear in the history view.



• ENDO MODE Horizontal and Vertical Periapicals


ENDO MODE can be used when you need multiple images of the same area during one appointment. While in ENDO MODE, single x-rays (either horizontal or vertical periapicals) are directly saved to the database. The image will go to the history view as well as the expanded view. To exit ENDO MODE, click on the stop button on the toolbar.

- ✓ BWX -- Four Horizontal
- BWX -- Four Vertical
- BWX -- Two Horizontal
- BWX -- Seven Vertical
- PA -- Single Horizontal
- PA -- Single Vertical
- Endo -- Horizontal
- Endo -- Vertical
- FMX -- Horizontal BWX
- FMX -- Eight Anterior
- FMX -- Vertical BWX
- Portrait -- Five Block




• Taking a layout series

Once you have selected the layout and **OK** button, the exposure screen will appear.




1. There will be a layout of boxes in the exposure screen.
 - 1.1. The numbers in the boxes indicate the order in which the x-rays will be exposed.
 - 1.2. If you wish to take the x-rays in a different order, you can use the arrow keys or the mouse to get you to the box where you want the x-ray to appear.
2. There will be a **red** outline around the box where the x-ray will be exposed.
3. The outline will change to **green** and you will hear a “ding” when the sensor is active and ready to be exposed.
4. Align X-Ray tube.
5. Push X-Ray button. You should hear another “ding” indicating that the exposure was detected.
6. Your image should appear in the active box within 2 - 6 seconds.
7. The active box will auto-advance to the next box in the series.
 - 7.1. If you wish to advance to a different box, just use the arrow keys or the mouse.
8. Repeat steps 4-7 until the layout is complete.
9. Save the images by single clicking on the **Save** icon on the tool bar.
10. After the images have been saved, they will appear in the history view.



• Reassigning Layouts to Patients

If a layout is exposed and saved under the wrong patient name, you can re-assign the layout by selecting it with “control+click” in the history view. In Lion (i.e. Mac OS 10.7.x) and later, a pop-over window will appear with the complete list of RV patients. Select the proper patient name and click the “Reassign” button to move the layout into the proper patient file.




• Deleting Patient Layouts

If a layout is no longer needed, it can be deleted from the database. To delete the layout, select it (or select a series using shift or command + click) in the history view, and press the “delete” key.

• Restoring an Exposure Series After a Failure to Save

When you choose to expose an xray layout or take an intra-oral camera image, RadioVision begins keeping track of all acquired images, and their locations in the selected layout. As these images are acquired, they are copied to a temporary location along with the selected layout and patient name. If for any reason RadioVision quits without saving this layout, the user will be presented with a dialog box allowing the user to recover this unsaved data.



If you choose to restore the data, the patient and layout will be reloaded. If this was an xray acquisition, the sensor will be re-armed for acquisition. If this was an IO camera series, the camera capture window will be opened. In either case, the user can continue with the data capture without data loss.

This is a very useful tool for recovering data that would otherwise be lost if a user chooses to select a new patient or quits RadioVision without saving the session.


If the previous session data is no longer needed, press the "Delete Last Session" to remove it. If the user begins a new acquisition, this session data will be deleted automatically, and replaced with the new session data.

• Retaking Images

Each location in a layout series or a single x-ray can hold multiple images stacked on top of each other.

If you need to retake an image, it must be done before the image or series of images has been saved.

1. Select the image that you need to retake with the arrow keys or the mouse.
2. The outline will change to **green** and you will hear a “ding” when the sensor is active and ready to be exposed.
3. Align X-Ray tube.
4. Push X-Ray button. You should hear another “ding” indicating that the exposure was detected.
5. Your image should appear in the active box within 6 seconds.
6. You will see (n/m) located in the upper left hand corner of the x-ray box. (n) is the nth image in a stack of (m) number of images.
7. To view other images in the stack hold down the option key and press the right or left arrow key. This will allow you to cycle through the images located in one box.
8. If you wish to delete an x-ray from the stack, simply press the delete key while the x-ray is highlighted.
9. Each stacked image can be opened separately in an expanded view and can have individual brightness and contrast settings applied.
10. Continue capturing the layout.
11. Save the images by single clicking on the **Save** icon on the tool bar.
12. After the images have been saved, they will appear in the history view.



• Hot-Swapping



If you have multiple sensor sizes, you may want to “hot-swap” (change sensors) in the middle of a layout.

1. When you remove the sensor from the USB box, there will be a **red** outline around the box where the x-ray will be exposed. This indicates that the sensor is not ready to be exposed.
2. Plug in the new sensor to the USB box.
3. The outline will change to **green** and you will hear a “ding” when the sensor is active and ready to be exposed.
4. Continue capturing the layout.


• Attaching notes to single images and layouts

Once you have saved images to the database you can attach notes to either the single image or to the entire layout.


1. When you run the cursor over the bottom of the layout in the history view, a button will appear. Single click on the **add notes** button. You can also run the cursor over the bottom of a single x-ray in the detail view and the **add notes** button will appear.




2. Single click on the [add notes](#) button. A window will appear and you will be able to add notes into the notes editor pane and view previous notes in the history pane. Once you click on the save button, the note will be put in the history pane and the date will be attached.



3. An image or layout with a note attached to it will be denoted in the lower left-hand corner. Click on the icon to view previous notes or to append notes.



Denotes note associated with layout



Denotes note associated with single x-ray

• Importing X-Rays

Once you are in a patient's chart, you are ready to begin importing x-rays

1. To import an x-ray, single click on the **Import** icon on the tool bar.
 - 1.1. Unlike taking an x-ray, no sensor needs to be plugged in while importing an image.
 - 1.2. A sheet will drop down from the tool bar and you will need to choose the type of x-rays you would like to import. Single click on the **layout type**
 - 1.3. A sheet will drop down you can select the layout you wish to import.
 - 1.3.1. If you are taking a single film, choose the correct orientation (vertical or horizontal). Then, select the teeth from the diagram that will be included in the film. You may select more than one tooth.
 - 1.3.2. If you are taking a layout series there is no need to select the teeth that will be included.
 - 1.3.3. If you have a series of images in one single image that you would like to import, choose the single horizontal PA. This will give you the largest image.
 - 1.4. You may edit the date to the date the image was taken.
 - 1.5. Once you have selected the layout, single click on the **OK** button.
 - 1.6. Simply drag and drop images into the appropriate box, or paste an image using command-v
 - 1.7. Save the image(s) by single clicking on the **Save** icon on the tool bar.

• Exporting X-Rays

Once you are in a patient's chart, you are ready to begin exporting x-rays

1. To export images to the desktop or directory, hold down the option key and drag the image from the detail view to your destination of choice.
2. Once the image is on the desktop it is in either TIFF or JPEG format and can be dragged and dropped into any graphics capable application (iPhoto, web based email etc)
3. To export an image to Apple's mail application, iWork or Text edit simply drag and drop the image from the detail view to the destination of choice. You can also use copy/paste (command-c/command-v) to copy from RV and paste into another application
4. If you would like to export a layout as a single image, drag the image from the history view.
5. Single images or layouts can also be copied onto the pasteboard and pasted into other applications

• Printing X-Rays

Once you are in a patient's chart, you are ready to begin printing x-rays


1. You can print images from the detail view or the expanded view.
2. If you would like to print the a layout, the images will need to be in the detail view
3. If you would like to print a single image from a layout, the image will need to be in the expanded view
4. Once you have selected the image(s) to be printed, select **File** from the menu bar.
5. A list will drop and you will select **Print**.
6. Another sheet will drop and you will need to select **Print** again.

To maximize the size of image printed you may want to select portrait or landscape from the page setup.

• Viewing X-Rays

Once you are in a patient's chart, you are ready to begin viewing x-rays

1. The history view contains all the x-rays taken on the patient in chronological order.
2. If the patient has had several x-rays taken you may need to use the scroll bar located at the bottom of the history view.
3. **DOUBLE CLICK** on the x-rays you would like to see and they will appear in the detail view.
4. If you would like to open multiple image layouts in the detail view, single click on the layouts you would like to open using the Shift and Command keys for multiple selections, and choose "Open Selected Layouts" from the file menu, or simply press command-O.





5. If you would like to enlarge or modify an x-ray image, **DOUBLE CLICK** on the x-ray in the detail view and it will appear in the expanded view.



6. To view image in full-screen mode, **DOUBLE CLICK** on the image in the expanded view.

• Filtering X-Rays by Tooth Number


1. Select “Filter by Tooth Number...” from the RadioVision File menu.
2. A sheet will drop down and you will select the teeth that you would like to view.
3. Click on the **OK** button.
4. X-rays of the teeth you selected will appear in the history view.
5. To restore the history view, select “Restore From Filtered” from the RadioVision File menu.



• Emailing X-Rays

1. Both the main RadioVision window and the expanded view window have an email toolbar item that can be used for pasting your radiographs into an email. RadioVision is configured to work with Apple's Mail application.
2. Use the main window to send images of entire layouts. RadioVision will create an image of the layout in the detail view, and pass it to the mail application.
3. Use the expanded view window to send emails of individual images within a layout

Clicking the email toolbar item in this configuration creates this email attachment



• Tools

RadioVision comes with a number of tools to enhance your images

1. In the DETAIL VIEW you are able to FLIP and ROTATE the images, and apply the CrystalView filter. All other tools can only be utilized in the expanded view.
2. Single click on the image that you would rotate, flip or CrystalView and click on the appropriate icon on the tool bar.
3. If you wish to save the new orientation of the image be sure to click on the [Save](#) icon in the tool bar after you have made your changes.
2. Adjust the BRIGHTNESS, CONTRAST and NOISE REDUCTION in the EXPANDED VIEW.
3. **Brightness** : uniformly changes the brightness of the image
4. **Contrast**: changes the image white and black points, so that bright spots become brighter, and dark spots become darker. This is useful for revealing subtle details in the image, where those details have very little change in brightness
5. **Noise Reduction**: moving the slider to the left increases the sharpness of lines in the image, and moving to the right softens lines.
6. When you are happy with the brightness, contrast and noise reduction be sure to click on the [Save](#) icon on the tool bar. These changes will be saved to the database.
7. You will always be able to revert to the original image by clicking on the [Revert to Original](#) icon on the toolbar.
8. **CrystalView Filter**: CrystalView is a Haas Software proprietary filter that improves a number of image characteristics simultaneously.



Original Image



Noise Reduction, Brightness
and Contrast Adjustments

- 8.1. CrystalView attempts to increase your image sharpness while also controlling the image noise. It also makes adjustments to the brightness, contrast and gamma correction. The settings for CrystalView are user-customizable, as described in the user preferences section of this manual.



Original Image




CrystalView Applied

9. **Highlight Filter** : The Highlight filter combines the effects of the CrystalView and False Color filters, helping to highlight features in a radiograph for the benefit of patient communication and education.



Original Image




Highlight Applied

10. **False Color Filter** : The false color filter applies the standard “rainbow mapping” to your grayscale radiograph. Dark areas are mapped to blue, bright areas to red, mid-point to green, and intensities throughout the range are given a combination of red, green and blue. The result is a much softer color image as compared to Highlight, but without the good definition and sharpening associated with Highlight.




Original Image




False Color Applied

11. **Backlight** : The Backlight filter produces some very interesting effects under certain conditions, by targeting bright objects in the radiograph. If a bright object contains no real information (i.e. the region is uniformly bright with no subtle variation) Backlight will turn the area quite dark. If a bright area does contain subtle changes in intensity that might indicate interest, that area is enhanced relative to the rest of the radiograph.




Original Image



Backlight Applied

12. **Equalize** : The Equalize filter performs what is known as a histogram equalization. In an unedited radiograph, the histogram distribution is typically tightly bunched, resulting in an image that does not have optimal contrast. Histogram equalization re-distributes the image intensities so that the histogram is completely flat: all intensities are equally represented. This filter is a bit harsh for many images, but if you acquire a particularly low-contrast image for some reason (or have one emailed to you), this filter will really pull out the details.



Original Image



Equalize Applied


13. **Magnify** : The magnify tool allows you to zoom in on different areas of an x-ray.

13.1. Click on the “magnify” toolbar item.

13.2. Click in the image area in the expanded x-ray view. This will cause the magnifying glass to appear, and will follow your cursor around the image. This is the “glass tracking mode”.

13.3. Click again in the image area. This will cause the magnifying glass to stick in place, and is now free from your cursor. This is the “glass viewing mode”.

13.4. Click in the image area again, but this time HOLD THE MOUSE BUTTON DOWN AND DRAGE THE MOUSE. This allows you to change the size and power of the magnifier. Dragging up and down changes the power; dragging left and right changes size. This is the “glass edit mode”. When you release the mouse button, you will again be in “glass viewing mode”. Clicking the mouse again will bring you to “glass viewing mode”.



Magnify Tool


14. **Measure** : The Measure tool allows you to measure the distance along a series of line segments.
- 14.1. Click on the “measure” toolbar item.
 - 14.2. Click on the image, at the top of the feature you’d like to measure. This sets the starting point of the first line segment.
 - 14.3. Move the mouse to the endpoint of the first line segment, and click. This creates a joint in the line defining the end of the first segment, and the beginning of the second.
 - 14.4. The length of each line segment is labelled, as is the total path length
 - 14.5. When you’ve finished the path, double-click the mouse to define the end of the last line segment, and terminate the path



Measure Tool

14.6. Calibrating your measurement

- 14.6.1. The scale of the measurement tool is initially set using the sensor pixel size. Primarily due to geometric foreshortening, the measure scale can sometimes be incorrect. RadioVision allows you to scale the measurement.
- 14.6.2. After you've terminated the measurement path, move the mouse to the bottom center of the image.
- 14.6.3. Click on the "calibrate measurement" button that appears.
- 14.6.4. Enter your calibration length (i.e. the known length of your path) in the measurement calibration panel, and click the "enter" button. Your known calibration length would typically come from a file inserted partially into the tooth.
- 14.6.5. If you click the "save" button on the expanded x-ray view toolbar, this calibration scale factor will be saved to the database, so that all future measurements will be calibrated and correct.



15. **Emboss** : The Emboss filter uses intensity differences in the image to create a 3D relief image. The level of emboss is adjustable, as described later in the user preferences section




Original Image



Embossed Image

• Intra-Oral Camera


1. RadioVision is compatible with multiple intra-oral cameras. Please contact your technical support for a current list of supported cameras.
2. The interface for all cameras is opened using the “open camera” toolbar item from the main window.
 - 2.1. When the window is opened, the camera begins acquiring images
 - 2.2. An image can be snapped using a number of methods
 - 2.2.1. Press the “snap” button in the image window
 - 2.2.2. Press the “start auto-snap” button in the image window, and hold the camera steady over the area of interest. RadioVision will automatically snap an image if you hold the camera still for more than a specified number of seconds, discussed later in the user preferences section.
 - 2.2.3. Support for the “snap” button on some cameras is supported. Contact your technical support for details.




• Creating Custom Layouts

RadioVision allows users to create custom layouts to suit the needs of the individual office. When a custom layout is created on one computer, that layout is available to all computers running RadioVision on the local network.

You start the process by choosing “manage” from the RadioVision Layouts menu item




Use “create new” to begin customizing your layout.



Clicking “Create New” in the Program Layouts window opens the Layout Creator window, shown below.


Use the row stepper to add or remove rows in each column.

Use the Column stepper to add or remove columns.




You start with one column, one row.

As you add columns and rows, RV will change the size and locations of the x-ray placeholders to optimize use of space.




Click on an individual place holder to change its orientation (horizontal or vertical exposure).




When the layout is complete, select the "Assign Teeth" button to disable the orientation changing mode, and enter tooth numbers for each location.

Select a location, then select the tooth numbers from the graphic below the layout.




The "Create" button will be enabled once all exposure locations have at least one tooth associated with them.

Name your layout before creating it.



• File Menu

Under the file menu you can backup your data, restore from the backup, enter your practice information and print.




When creating a backup with the “Backup Now” File menu command, the entire contents of the database are copied to a file in a location that you specify. When restoring from a backup using the “Restore From Backup...” File menu command, all users must be disconnected from the database. When creating backups, or restoring from a backup using the aforementioned File menu commands, PostgreSQL requires that the user logged in must also be the name of a PostgreSQL database user.

• Remote Data Backup

Users with active maintenance agreements have access to the Haas Software remote backup server. The remote backup provides two major advantages over the local database backup.

1. Your data is continuously and incrementally backed up. Within 60 seconds of saving new data, the remote backup feature copies this new data to the remote server so that you are completely protected against hardware failures.
2. You have fine-grain control over the data you want to recover. Consider the following as an example: Your local backup is set to occur at 6 pm. At 4 pm, someone accidentally deletes a full-mouth series from a patient record. If you only have the local RadioVision backup, you will have to replace the entire contents of your database in order to recover the deleted layout; this will result in losing all layouts that were taken since the

last local backup, which would have been 6 pm the previous evening. If you have data backed up to the remote server, you can simply select the deleted layout from the list of available layouts and recover it; no data will be lost.




Selecting “Remote Backup Status..” from the File menu displays a summary of your protected data in the “RadioVision Remote Storage” window. The “Remote

The screenshot shows the RadioVision Remote Storage window. It contains a table with the following data:


Table Name	Backup Up	Not Backup Up	Table Status
Patients	2	0	Idle
Layouts	8	0	Idle
Layout Templates	1	0	Idle
Images	32	0	Idle
Practice Info	1	0	Idle

Remote Server Status: Connected


Server Status” text on the lower left portion of the “Remote Storage Status” window will indicate whether the remote server is connected. If you are connected but do not have an active maintenance agreement, the status text will indicate this.




Select “Recover Missing Data...” from the File menu to recover data from the remote server, which opens the “Data Recovery” window.




The tables in the Data Recovery window display RECOVERABLE DATA ONLY! If all of the data on the remote backup server also exists in your local database (i.e. you haven't deleted any data since creating it), there will be NO recoverable data. If you delete data from your local database, it will be displayed as recoverable data from the remote backup.



In order to maintain a unique relationship between all data objects, a strict one-to-one relationship between the local and remote databases is enforced. If you restore a previous local database backup using the "Restore From Backup..." File menu command, a new remote database will be created using the creation date information contained in your local database. The "Remote Databases" pop-up button shows the list of databases available. The remote database that corresponds to your current local database will have the parenthetic phrase "Merge With" displayed, indicating that your recovered data will be merged with your local database. Remote databases that do not correspond to your local database will have the parenthetic phrase "Replace With" displayed, indicating that your recovered data will replace data in your database, and your local database will from that time forward be linked to this chosen remote database.




The “Data Recovery” window above shows the result of deleting the Demo Patient from the database, and then displaying the recoverable data. Deleting a patient from the database also deletes layouts and images associated with that patient. In order to recover that data, you must select the data in each list, and then press the “Recover Selected Data” button. If you only select “Demo Patient” from the recoverable patients table, but do not select any layouts or images, you will not recover that patient’s image data.



• Preferences

RadioVision Preferences allow each user with separate Mac OS X user account to customize their settings.



1. General Tab

1.1. Print practice information on images


- 1.1.1. Causes text you enter in into the practice information to be printed on all RadioVision printouts

1.2. Size of Median Filter


- 1.2.1. Sets the size of the nxn median filter box. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED
- 1.2.2. A median filter setting of 1 disables the median filter

1.3. Emboss Level

- 1.3.1. Sets the degree of emboss. A higher number is a more dramatic effect



- 1.4. Data Smoothing Level
 - 1.4.1. An alternate approach to median filtering. Data smoothing is a more sophisticated algorithm than the median filter, accounting for image gradients within the nxn box. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED.
 - 1.4.2. A smoothing level of 1 disables the data smoothing filter
 - 1.5. Apply CrystalView to Acquired Images
 - 1.5.1. Causes the CrystalView filter to be applied to acquired images automatically. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED. This filter CANNOT BE UNDONE, AND PERMANENTLY ALTERS THE IMAGE DATA.
 - 1.6. Apply Highlight to Acquired Images
 - 1.6.1. Causes the Highlight filter to be applied to acquired images automatically. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED. This filter CANNOT BE UNDONE, AND PERMANENTLY ALTERS THE IMAGE DATA.
 - 1.7. Use TIFF/JPEG for exported x-rays or database storage
 - 1.7.1. Causes RadioVision to export or save database images in either TIFF or JPEG format. The TIFF format is uncompressed, and there is no loss of data, but the files are quite large. JPEG is a lossy compression technique, but does produce high-quality, smaller images.
 - 1.8. Play Sound on Acquisition
 - 1.8.1. A “ping” will sound when the x-ray sensor has detected an exposure, if selected.
 - 1.9. Play Sound on Ready to Acquire
 - 1.9.1. A “ding” will sound when the x-ray sensor is ready to be exposed, if selected
2. Backup Tab
 - 2.1. Auto Backup
 - 2.1.1. If selected, RadioVision will create a database backup at the time indicated by “Daily Backup Time”
 - 2.2. Path to Backup File
 - 2.2.1. Indicates the folder in which the backup file will be created. You may enter a path manually, or use the “Browse...” button to navigate to the desired path.
 - 2.3. Number of Backups to Keep
 - 2.3.1.1. RadioVision will keep the last “n” backup files. This gives you the option to keep a history of backups, while keeping the total backup size small. The backup file will be



named RadioVision_YYYY-mm-dd_hh_MM.back, where y=year, m=month, d=day, h=hour, M=minute.

2.4. Compress backup files

2.4.1. If selected, the backup files are compressed after being created which can be useful if these backup files are being backed up by Time Machine or an equivalent.

2.5. Enable Remote Backup

2.5.1. If selected, RadioVision will connect to remote backup server, and automatically back up all data.


2.6. Remote Backup Server Address

2.6.1. The address used to connect to the remote backup server.

3. Layout Options Tab


3.1. Pop-up menu on the left allows you to select the layout you'd like to modify

3.1.1. When you choose a layout, a series of layout locations labeled with numbers appears. The numbers indicate the exposure sequence for that layout. To edit the the layout, do the following:




Note: When taking a layout series you can always deviate from the automated order by highlighting the frame that you want to expose.


- 3.1.1.1. Press the “Clear” button
- 3.1.1.2. Click the layout locations in the order you’d like to expose them
 - 3.1.1.2.1. Click a location a second time in order to remove the sequence number, in the event that the wrong location is accidentally clicked.
- 3.1.1.3. Click the “Apply” button when the sequence is complete.
- 3.1.2. Click the “Restore Defaults” location to re-apply the default sequence for the displayed layout.
- 3.2. Pop-up menu on the right allows you to specify which layout should be selected by default when acquiring or importing x-rays.
4. Server Tab
 - 4.1. Use Bonjour, or Use Fixed IP Address




- 4.1.1. As the default RV uses Bonjour to establish client connections to the server. Improper network configuration can cause this automated process to fail. As an alternative to Bonjour you can choose to use a fixed server IP Address which avoids the need for domain name server lookups.
- 4.2. Automatically log in as user:
 - 4.2.1. In order to avoid the need to type in the user name and password each time RadioVision starts up, you can select this option and enter the information in this panel.
IMPORTANT: YOU MUST HIT THE "RETURN" KEY AFTER ENTERING THE PASSWORD IN ORDER FOR THE PASSWORD INFORMATION TO BE PASSED TO THE USER DEFAULTS!
- 5. Sensor Tab
 - 5.1. Sensitivity
 - 5.1.1. A higher number indicates a higher sensitivity setting for the sensor. A higher sensitivity, or gain, is at the expense of higher noise.
 - 5.1.1.1. This option is only enabled when using Suni sensors.
 - 5.2. Integration Time (ms)
 - 5.2.1. Sets the sensor integration time (time during which charge created from x-ray exposure is collected by the sensor).
 - 5.2.1.1. This option is only enabled when using Suni sensors.
 - 5.3. Threshold
 - 5.3.1. Sets the minimum charge threshold used to trigger an exposure
 - 5.3.1.1. This option is only enabled when using Suni sensors.
 - 5.4. Flush delay (s)
 - 5.4.1. Sets the number of seconds the sensor will spend flushing charge before going to the "ready to expose" state.
 - 5.4.1.1. This option is only enabled when using Suni sensors.
 - 5.5. Stretch



- 5.5.1. This is a non-linear contrast enhancement setting used to improve the image histogram. The image intensities at the higher and lower ends are re-distributed to more fully populate the image histogram.
- 5.5.2. A lower stretch number is a more aggressive re-distribution, where a setting of 5 essentially disables this feature. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED.
- 5.6. Pedestal
 - 5.6.1. This determines the level of image noise rejection. A higher pedestal numbers offers better noise rejection, and improves image contrast by ensuring the image histogram is populated by real signal, not noise. THIS SETTING IS APPLIED TO IMAGES AS THEY ARE ACQUIRED, AND HAS NO EFFECT ON IMAGES ALREADY ACQUIRED.
- 6. Filters Tab
 - 6.1. Pop-up menu allows you to choose the filter you'd like to customize
 - 6.2. CrystalView settings
 - 6.2.1. Radius: the local extent to which the filter is applied. The filter is applied to the entire image, but a larger radius means more pixels are used to determine the filtered value of any particular pixel
 - 6.2.2. Sharpness: sets how aggressively RadioVision sharpens the image. As with all sharpening filters, more sharpness can lead to more noise or graininess



- 6.2.3. Level: sets the mid-tone brightness of the filtered image
- 6.3. Additional Highlight settings
 - 6.3.1. Smooth: sets the smoothing to be applied (noise reduction) to the image prior to applying the remaining filter values
- 6.4. Restore Defaults
 - 6.4.1. Re-applies the RadioVision default values for the selected filter
- 6.5. Upper right image pane
 - 6.5.1. Drop an un-filtered image into this pane to see how your filter settings will behave
- 6.6. Lower right image pane (with the RadioVision logo by default)
 - 6.6.1. Your filter settings are applied in real-time, and this lower pane shows the filtered version of the image you dropped into the upper frame.
- 7. Advanced Tab
 - 7.1. Suni Plus Vendor ID
 - 7.1.1. If your Suni Plus sensor is not recognized when you plug it in, there is a chance that a bit in the EEPROM indicating the vendor ID has been flipped from 1 to 0. Changing the expected ID accordingly will allow RadioVision to once again recognize your hardware. **DO NOT CHANGE THIS SETTING WITHOUT THE DIRECTION OF YOUR TECHNICAL SUPPORT**
 - 7.2. Camera Movement Threshold
 - 7.2.1. When using the intra-oral camera, RadioVision will “auto-snap” an image if you dwell over one area long enough. The radius within which RadioVision considers that the image has not moved is set by this threshold. If you’re having a hard time holding the camera still enough to capture an image with “auto-snap”, try increasing this threshold.
 - 7.3. Camera Time to Snap
 - 7.3.1. Sets the time threshold during which the intra-oral camera image must remain within the “camera movement threshold” radius before an image will be auto-snapped.
 - 7.4. Enable Demo Sensor
 - 7.4.1. If selected, you will have the option of using a software sensor to acquire x-rays. This is useful if you want to see how the acquisition process works, but don’t have x-ray



hardware available. If selected, the DemoSize2 sensor will be available in the “sensor” list in the sheet that appears when you choose to take an x-ray.

7.5. Enable QuickRay Sensor


7.5.1. Please contact your technical support if you are interested in QuickRay sensor support

7.6. License Server Port

7.6.1. This controls the port number used to communicate license information from the RadioVision server to the clients. This should normally be set to “0”, which allows RadioVision to find an open port number, and communicate that information to the clients using the Bonjour service discovery. If you are using clients that are not on the same LAN as your server (i.e. you have a server in one office, and clients in a different office) you must set this server port to a non-zero number and use port mapping on your server LAN router to forward that port traffic to your server.


• Troubleshooting

1. RadioVision starts up, but the **select patient** button is disabled, and it is not searching for the server.



1.1. If this is a client computer




- 1.1.1. Verify that RadioVision is running on the server.
- 1.1.2. If the server is not on and running RadioVision, quit RadioVision on the client, turn the server on and start RadioVision on the server, then start RadioVision on the client.
- 1.1.3. If the problem persists, select **Re-Find Server** from the pop-up list in the lower left hand corner of the application window.



1.2. If this is a server computer

1.2.1. Select **Re-Find Server** from the pop-up list in the lower left hand corner of the application window.

1.2.2. The program will begin searching for the server and you will need to click on the **Use Local** button.



2. Image quality questions

2.1. For questions relating to sensor performance, please contact your technical support representative.